

C17orf64 Antibody, HRP conjugated



PACO37371

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	C17orf64
Source:	Uniprot
Rabbit	Q86WR6
Isotype:	Synonyms:
IgG	Uncharacterized protein C17orf64, C17orf64
Applications:	Immunogen:
ELISA	Recombinant Human Uncharacterized protein C17orf64 protein (111-236AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A