

ABI3BP Antibody, HRP conjugated



PACO37395

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular space.

Gene ID:

ABI3BP

Uniprot

Q7Z7G0

Synonyms:

Target of Nesh-SH3 (Tarsh) (ABI gene family member 3-binding protein) (Nesh-binding protein) (NeshBP), ABI3BP, NESHBP TARSH

Immunogen:

Recombinant Human Target of Nesh-SH3 protein (601-900AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A