

LYG2 Antibody



PACO37402

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May act as a potent antibacterial protein that may play a role in the innate immunity.

Gene ID:

LYG2

Uniprot

Q86SG7

Synonyms:

Lysozyme g-like protein 2 (EC 3.2.1. -), LYG2, LYGH


Immunogen:

Recombinant Human Lysozyme g-like protein 2 protein (20-212AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO37402 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37402 at dilution of 1:100.