

GPR141 Antibody, FITC conjugated



PACO37452

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Orphan receptor.

Gene ID:

GPR141

Uniprot

Q7Z602

Synonyms:

Probable G-protein coupled receptor 141 (G-protein coupled receptor PGR13), GPR141, PGR13

Immunogen:

Recombinant Human Probable G-protein coupled receptor 141 protein (153-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A