

PDZD3 Antibody, FITC conjugated



PACO37468

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a regulatory protein that associates with GUCY2C and negatively modulates its heat-stable enterotoxin-mediated activation. Stimulates SLC9A3 activity in the presence of elevated calcium ions.

Gene ID:

PDZD3

Uniprot

Q86UT5

Synonyms:

Na(+)/H(+) exchange regulatory cofactor NHE-RF4 (NHERF-4) (Intestinal and kidney-enriched PDZ protein) (Natrium-phosphate cotransporter IIa C-terminal-associated protein 2) (Na/Pi cotransporter C-terminal-associated protein 2) (NaPi-Cap2) (PDZ domain-containing protein 2) (PDZ domain-containing protein 3) (Sodium-hydrogen exchanger regulatory factor 4), PDZD3, IKEPP NHERF4 PDZK2

Immunogen:

Recombinant Human Na(+)/H(+) exchange regulatory cofactor NHE-RF4 protein (42-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A