

# FAM131B Antibody



PACO37478

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FAM131B

**Uniprot**

Q86XD5

**Synonyms:**

Protein FAM131B, FAM131B, KIAA0773

**Immunogen:**

Recombinant Human Protein FAM131B protein (1-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A