

# VMO1 Antibody



PACO37502

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

**Gene ID:**

VMO1

**Uniprot**

Q7Z5L0

**Synonyms:**

Vitelline membrane outer layer protein 1 homolog, VMO1

**Immunogen:**


Recombinant Human Vitelline membrane outer layer protein 1 homolog protein (25-202AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images


---



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO37502 at dilution of 1:100.



Immunofluorescent analysis of Hela cells using PACO37502 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human liver tissue using PACO37502 at dilution of 1:100.