

PACO37510

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Type I interferon required for maintaining basal levels of IFN-regulated genes, including 2'-5'-oligoadenylate synthetase, IRF7 and ISG15, in the female reproductive tract. Directly mediates protection against viral, including HSV-2, and bacterial, including Chlamydia muridarum, genital infections.

**Gene ID:**

Ifne

**Uniprot**

Q80ZF2

**Synonyms:**

Interferon epsilon (IFN-epsilon) (Interferon epsilon-1) (Interferon tau-1), Ifne, Ifne1 Ifnt1

**Immunogen:**

Recombinant Mouse Interferon epsilon protein (22-192AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A