

Ifne Antibody, HRP conjugated



PACO37511

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Type I interferon required for maintaining basal levels of IFN-regulated genes, including 2'-5'-oligoadenylate synthetase, IRF7 and ISG15, in the female reproductive tract. Directly mediates protection against viral, including HSV-2, and bacterial, including Chlamydia muridarum, genital infections.

Gene ID:

Ifne

Uniprot

Q80ZF2

Synonyms:

Interferon epsilon (IFN-epsilon) (Interferon epsilon-1) (Interferon tau-1), Ifne, Ifne1 Ifnt1

Immunogen:

Recombinant Mouse Interferon epsilon protein (22-192AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A