

FAM124A Antibody



PACO37518

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FAM124A

Uniprot

Q86V42

Synonyms:

Protein FAM124A, FAM124A


Immunogen:

Recombinant Human Protein FAM124A protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO37518 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO37518 at dilution of 1:100.