

# LDHD Antibody, Biotin conjugated



PACO37529

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

mitochondrion.

**Gene ID:**

LDHD

**Uniprot**

Q86WU2

**Synonyms:**

Probable D-lactate dehydrogenase, mitochondrial (DLD) (Lactate dehydrogenase D) (EC 1.1.2.4), LDHD

**Immunogen:**

Recombinant Human Probable D-lactate dehydrogenase, mitochondrial protein (310-507AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A