

MAMDC2 Antibody



PACO37534

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MAMDC2

Uniprot

Q7Z304

Synonyms:

MAM domain-containing protein 2 (MAM domain-containing proteoglycan) (Mamcan), MAMDC2

Immunogen:

Recombinant Human MAM domain-containing protein 2 protein (455-686AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A