

VN1R5 Antibody



PACO37546

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Putative pheromone receptor.

Gene ID:

VN1R5

Uniprot

Q7Z5H4

Synonyms:

Vomer nasal type-1 receptor 5 (G-protein coupled receptor GPCR26) (hGPCR26) (V1r-like receptor 5), VN1R5, V1RL5


Immunogen:

Recombinant Human Vomer nasal type-1 receptor 5 protein (185-269AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO37546 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO37546 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).