

TRAPPC6B Antibody, HRP conjugated



PACO37559

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in vesicular transport from endoplasmic reticulum to Golgi.

Gene ID:

TRAPPC6B

Uniprot

Q86SZ2

Synonyms:

Trafficking protein particle complex subunit 6B, TRAPPC6B

Immunogen:

Recombinant Human Trafficking protein particle complex subunit 6B protein (1-130AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A