

# RTN4RL1 Antibody, FITC conjugated



PACO37564

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

**Gene ID:**

RTN4RL1

**Uniprot**

Q86UN2

**Synonyms:**

Reticulon-4 receptor-like 1 (Nogo receptor-like 2) (Nogo-66 receptor homolog 2) (Nogo-66 receptor-related protein 3) (Ngr3), RTN4RL1, NGRH2 NGRL2

**Immunogen:**

Recombinant Human Reticulon-4 receptor-like 1 protein (201-420AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A