

RTN4RL1 Antibody, Biotin conjugated



PACO37565

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

Gene ID:

RTN4RL1

Uniprot

Q86UN2

Synonyms:

Reticulon-4 receptor-like 1 (Nogo receptor-like 2) (Nogo-66 receptor homolog 2) (Nogo-66 receptor-related protein 3) (Ngr3), RTN4RL1, NGRH2 NGRL2

Immunogen:

Recombinant Human Reticulon-4 receptor-like 1 protein (201-420AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A