

ZSCAN29 Antibody



PACO37614

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZSCAN29

Uniprot

Q8IWY8

Synonyms:

Zinc finger and SCAN domain-containing protein 29 (Zinc finger protein 690), ZSCAN29, ZNF690


Immunogen:

Recombinant Human Zinc finger and SCAN domain-containing protein 29 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human melanoma using PACO37614 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO37614 at dilution of 1:100.