

C11orf53 Antibody



PACO37618

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

C11orf53

Uniprot

Q8IXP5

Synonyms:

Uncharacterized protein C11orf53, C11orf53


Immunogen:

Recombinant Human Uncharacterized protein C11orf53 protein (13-236AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: C11orf53 antibody at 0.8 μ g/ml + Mouse gonadal tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 26 kDa. Observed band size: 26 kDa.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO37618 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37618 at dilution of 1:100.