

SEPT12 Antibody, Biotin conjugated



PACO37633

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Filament-forming cytoskeletal GTPase. May play a role in cytokinesis (Potential).

Gene ID:

41153

Uniprot

Q8IYM1

Synonyms:

Septin-12, SEPT12

Immunogen:

Recombinant Human Septin-12 protein (1-358AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A