

UBE2DNL Antibody



PACO37666

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

UBE2DNL

Uniprot

Q8IWF7

Synonyms:

Putative ubiquitin-conjugating enzyme E2 D2-like protein (Ubiquitin carrier protein D2-like) (Ubiquitin-conjugating enzyme E2D N-terminal-like) (Ubiquitin-protein ligase D2-like), UBE2DNL, UBE2D2L

Immunogen:

Recombinant Human Putative ubiquitin-conjugating enzyme E2 D2-like protein (1-75AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A