

RTKN2 Antibody, HRP conjugated



PACO37683

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an important role in lymphopoiesis.

Gene ID:

RTKN2

Uniprot

Q8IZC4

Synonyms:

Rhotekin-2 (Pleckstrin homology domain-containing family K member 1) (PH domain-containing family K member 1), RTKN2, PLEKHK1

Immunogen:

Recombinant Human Rhotekin-2 protein (1-163AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A