

NUDCD3 Antibody, Biotin conjugated



PACO37697

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NUDCD3

Uniprot

Q8IVD9

Synonyms:

NudC domain-containing protein 3, NUDCD3, KIAA1068

Immunogen:

Recombinant Human NudC domain-containing protein 3 protein (142-361AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A