

UBQLNL Antibody



PACO37718

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

UBQLNL

Uniprot

Q8IYU4

Synonyms:

Ubiquilin-like protein, UBQLNL


Immunogen:

Recombinant Human Ubiquilin-like protein (305-475AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37718 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37718 at dilution of 1:100.