

ABCA9 Antibody



PACO37730

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May play a role in monocyte differentiation and lipid homeostasis.

Gene ID:

ABCA9

Uniprot

Q8IUU7

Synonyms:

ATP-binding cassette sub-family A member 9, ABCA9


Immunogen:

Recombinant Human ATP-binding cassette sub-family A member 9 protein (106-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO37730 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).