

C20orf196 Antibody



PACO37742

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C20orf196

Uniprot

Q8IYI0

Synonyms:

Uncharacterized protein C20orf196, C20orf196

Immunogen:

Recombinant Human Uncharacterized protein C20orf196 protein (1-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A