

ODF3L1 Antibody



PACO37786

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

ODF3L1

Uniprot

Q8IXM7

Synonyms:

Outer dense fiber protein 3-like protein 1, ODF3L1


Immunogen:

Recombinant Human Outer dense fiber protein 3-like protein 1 protein (1-274AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: ODF3L1 antibody at 3 μ g/ml + HT29 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 32 kDa. Observed band size: 32 kDa.



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO37786 at dilution of 1:100.