

# KRTAP13-1 Antibody



PACO37806

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins.

**Gene ID:**

KRTAP13-1

**Uniprot**

Q8IUUC0

**Synonyms:**

Keratin-associated protein 13-1 (High sulfur keratin-associated protein 13.1), KRTAP13-1, KAP13.1 KRTAP13.1

**Immunogen:**

Recombinant Human Keratin-associated protein 13-1 protein (1-172AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A