

SLITRK4 Antibody, Biotin conjugated



PACO37817

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Suppresses neurite outgrowth.

Gene ID:

SLITRK4

Uniprot

Q8IW52

Synonyms:

SLIT and NTRK-like protein 4, SLITRK4

Immunogen:

Recombinant Human SLIT and NTRK-like protein 4 protein (201-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A