

# PSTK Antibody, HRP conjugated



PACO37831

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Specifically phosphorylates seryl-tRNA(Sec) to O-phosphoseryl-tRNA(Sec), an activated intermediate for selenocysteine biosynthesis.

**Gene ID:**

PSTK

**Uniprot**

Q8IV42

**Synonyms:**

L-seryl-tRNA(Sec) kinase (EC 2.7.1.164) (O-phosphoseryl-tRNA(Sec) kinase), PSTK, C10orf89

**Immunogen:**

Recombinant Human L-seryl-tRNA(Sec) kinase protein (1-347AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A