

PSTK Antibody, FITC conjugated



PACO37832

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Specifically phosphorylates seryl-tRNA(Sec) to O-phosphoseryl-tRNA(Sec), an activated intermediate for selenocysteine biosynthesis.

Gene ID:

PSTK

Uniprot

Q8IV42

Synonyms:

L-seryl-tRNA(Sec) kinase (EC 2.7.1.164) (O-phosphoseryl-tRNA(Sec) kinase), PSTK, C10orf89

Immunogen:

Recombinant Human L-seryl-tRNA(Sec) kinase protein (1-347AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A