
Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Multivalent adapter protein which may decrease NOS3 activity by inducing its translocation away from the plasma membrane.

Gene ID:

NOSTRIN

Uniprot

Q8IVI9

Synonyms:

Nostrin (BM247 homolog) (Nitric oxide synthase traffic inducer) (Nitric oxide synthase trafficker) (eNOS-trafficking inducer), NOSTRIN


Immunogen:

Recombinant Human Nostrin protein (1-322AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: NOSTRIN antibody at 6 μ g/ml. Lane 1: Mouse heart tissue. Lane 2: Mouse brain tissue. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 58, 55, 49, 65 kDa. Observed band size: 58 kDa.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO37854 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO37854 at dilution of 1:100.