

TIGD4 Antibody



PACO37882

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

TIGD4

Uniprot

Q8IY51

Synonyms:

Tigger transposable element-derived protein 4, TIGD4


Immunogen:

Recombinant Human Tigger transposable element-derived protein 4 protein (301-512AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37882 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skin tissue using PACO37882 at dilution of 1:100.