

# TAGAP Antibody



PACO37910

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May function as a GTPase-activating protein and may play important roles during T-cell activation.

**Gene ID:**

TAGAP

**Uniprot**

Q8N103

**Synonyms:**

T-cell activation Rho GTPase-activating protein (T-cell activation GTPase-activating protein), TAGAP, TAGAP1

**Immunogen:**


Recombinant Human T-cell activation Rho GTPase-activating protein (97-260AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO37910 at dilution of 1:100.