

# C10orf53 Antibody, HRP conjugated



PACO37947

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C10orf53

**Uniprot**

Q8N6V4

**Synonyms:**

UPF0728 protein C10orf53, C10orf53

**Immunogen:**

Recombinant Human UPF0728 protein C10orf53 protein (46-116AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A