

# C11orf45 Antibody, Biotin conjugated



PACO38009

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C11orf45

**Uniprot**

Q8TAV5

**Synonyms:**

Putative uncharacterized protein C11orf45, C11orf45

**Immunogen:**

Recombinant Human Putative uncharacterized protein C11orf45 protein (23-145AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A