

DNAJB14 Antibody, HRP conjugated



PACO38015

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a co-chaperone with HSPA8 and promotes the degradation of misfolded transmembrane proteins in the ER-associated degradation (ERAD) pathway.

Gene ID:

DNAJB14

Uniprot

Q8TBM8

Synonyms:

DnaJ homolog subfamily B member 14, DNAJB14

Immunogen:

Recombinant Human DnaJ homolog subfamily B member 14 protein (1-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A