

# MRPL30 Antibody



PACO38022

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

mitochondrial inner membrane, mitochondrion, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

**Gene ID:**

MRPL30

**Uniprot**

Q8TCC3

**Synonyms:**

39S ribosomal protein L30, mitochondrial (L30mt) (MRP-L30) (39S ribosomal protein L28, mitochondrial) (L28mt) (MRP-L28) (Mitochondrial large ribosomal subunit protein uL30m), MRPL30, MRPL28 RPML28

**Immunogen:**

Recombinant Human 39S ribosomal protein L30, mitochondrial protein (35-161AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A