

PACO38062

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.

**Gene ID:**

LRRC29

**Uniprot**

Q8WV35

**Synonyms:**

Leucine-rich repeat-containing protein 29 (F-box and leucine-rich repeat protein 9) (F-box protein FBL9) (F-box/LRR-repeat protein 9), LRRC29, FBL9 FBXL9

**Immunogen:**

Recombinant Human Leucine-rich repeat-containing protein 29 protein (1-223AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A