

ST8SIA4 Antibody, HRP conjugated



PACO38083

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the polycondensation of alpha-2,8-linked sialic acid, required for the synthesis of polysialic acid, (PSA), which is present on the embryonic neural cell adhesion molecule (N-CAM), necessary for plasticity of neural cells.

Gene ID:

ST8SIA4

Uniprot

Q92187

Synonyms:

CMP-N-acetylneuraminate-poly-alpha-2,8-sialyltransferase (EC 2.4.99. -) (Alpha-2,8-sialyltransferase 8D) (Polysialyltransferase-1) (Sialyltransferase 8D) (SIAT8-D) (Sialyltransferase St8Sia IV) (ST8SiaIV), ST8SIA4, PST PST1 SIAT8D

Immunogen:

Recombinant Human CMP-N-acetylneuraminate-poly-alpha-2,8-sialyltransferase protein (21-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A