

ATG2B Antibody, Biotin conjugated



PACO38121

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for both autophagosome formation and regulation of lipid droplet morphology and dispersion.

Gene ID:

ATG2B

Uniprot

Q96BY7

Synonyms:

Autophagy-related protein 2 homolog B, ATG2B, C14orf103

Immunogen:

Recombinant Human Autophagy-related protein 2 homolog B protein (1357-1536AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A