

CXCL10 Antibody, Biotin conjugated



PACO38209

Product Information

Size:

50ug

Reactivity:

Macaca mulatta

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Chemotactic for monocytes and T-lymphocytes. Binds to CXCR3.

Gene ID:

CXCL10

Uniprot

Q8MIZ1

Synonyms:

C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein) (Gamma-IP10) (IP-10) (Small-inducible cytokine B10), CXCL10, SCYB10

Immunogen:

Recombinant Macaca mulatta C-X-C motif chemokine 10 protein (22-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A