

C20orf197 Antibody



PACO38230

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C20orf197

Uniprot

Q8N268

Synonyms:

Uncharacterized protein C20orf197, C20orf197

Immunogen:

Recombinant Human Uncharacterized protein C20orf197 protein (1-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A