

STX19 Antibody



PACO38254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Plays a role in endosomal trafficking of the epidermal growth factor receptor (EGFR).

Gene ID:

STX19

Uniprot

Q8N4C7

Synonyms:

Syntaxin-19, STX19


Immunogen:

Recombinant Human Syntaxin-19 protein (1-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: STX19 antibody at 2 μ g/ml + K562 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 35 kDa. Observed band size: 35 kDa.