

SPEM1 Antibody, Biotin conjugated



PACO38261

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for proper cytoplasm removal during spermatogenesis.

Gene ID:

SPEM1

Uniprot

Q8N4L4

Synonyms:

Spermatid maturation protein 1, SPEM1, C17orf83

Immunogen:

Recombinant Human Spermatid maturation protein 1 protein (50-309AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A