

C12orf50 Antibody, FITC conjugated



PACO38304

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

C12orf50

Uniprot

Q8NA57

Synonyms:

Uncharacterized protein C12orf50, C12orf50

Immunogen:

Recombinant Human Uncharacterized protein C12orf50 protein (201-414AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A