

SIMC1 Antibody



PACO38318

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

SUMO polymer binding.

Gene ID:

SIMC1

Uniprot

Q8NDZ2

Synonyms:

SUMO-interacting motif-containing protein 1, SIMC1, C5orf25

Immunogen:

Recombinant Human SUMO-interacting motif-containing protein 1 protein (18-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A