

NSUN7 Antibody



PACO38322

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May have S-adenosyl-L-methionine-dependent methyl-transferase activity.

Gene ID:

NSUN7

Uniprot

Q8NE18

Synonyms:

Putative methyltransferase NSUN7 (EC 2.1.1. -) (NOL1/NOP2/Sun domain family member 7), NSUN7


Immunogen:

Recombinant Human Putative methyltransferase NSUN7 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO38322 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO38322 at dilution of 1:100.