

# LSM10 Antibody, Biotin conjugated



PACO38405

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	LSM10
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q969L4
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	U7 snRNA-associated Sm-like protein LSm10, LSM10
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human U7 snRNA-associated Sm-like protein LSm10 protein (1-123AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A