

TRIM43 Antibody



PACO38422

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:**

TRIM43

Uniprot

Q96BQ3

Synonyms:

Tripartite motif-containing protein 43, TRIM43


Immunogen:

Recombinant Human Tripartite motif-containing protein 43 protein (43-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot. Positive WB detected in: Mouse brain tissue. All lanes: TRIM43 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 53 kDa. Observed band size: 53 kDa.