

TTC9C Antibody



PACO38546

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

cytoplasm, endoplasmic reticulum membrane, FK506 binding, peptidyl-prolyl cis-trans isomerase activity, chaperone-mediated protein folding, protein peptidyl-prolyl isomerization.

Gene ID:

TTC9C

Uniprot

Q8N5M4

Synonyms:

Tetratricopeptide repeat protein 9C (TPR repeat protein 9C), TTC9C


Immunogen:

Recombinant Human Tetratricopeptide repeat protein 9C protein (1-171AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: TTC9C antibody at 12 μ g/ml + THP-1 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 21, 25 kDa. Observed band size: 21 kDa.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO38546 at dilution of 1:100.